

Controlled temperature calibration micro baths

LR-Cal FLUID 100-N: -18...+150°C

LR-Cal FLUID 100-45: -30...+150°C

(stated temperature range valid at ambient temperature+20°C)

The portable temperature calibration micro baths **LR-Cal FLUID 100-N** and **LR-Cal FLUID 100-45** serve as temperature source and reference instrument in one. For testing, adjusting and calibrating all types of temperature measuring instruments. They are also particularly suitable for laboratory and glass thermometer as well as temperature probes with e.g. 90° bends.



Technical data:

Model:	LR-Cal FLUID 100-N	LR-Cal FLUID 100-45
Temperature range at +20°C ambient temperature:	-18...+150°C	-30...+150°C
Accuracy of the temperature indication:	±0.15°C	±0.015°C
Display resolution (switchable):	0.01° / 0.1°	0.01° / 0.1°
Temperature units (selectable):	°C / °F / K	°C / °F / K
Stability of regulated temperature:	±0.02°C (1)	±0.02°C (1)
Average heating time:	12°C per minute	3.5°C per minute
Average cooling time:	5°C per minute	1°C per minute
Radial temperature uniformity:	±0.04°C (2)	±0.04°C (2)
Axial temperature uniformity (in the holes):	±0.05°C	±0.05°C
Reservoir depth:	185 mm (usable: 150 mm)	185 mm (usable: 150 mm)
Reservoir diameter:	54 mm	45 mm
Display:	LED, 2 lines	LED, 2 lines
Interface:	RS232 (Option: USB-converter)	RS232 (Option: USB-converter)
Inputs for 2 external probes: (3)	Option (version -2I)	Option (version -2I)
Ramp function (slope):	•	•
Thermostat test function and connection:	•	•
Power supply (50/60 Hz):	230 VAC (4)	230 VAC
Power consumption:	300 VA	300 VA
Weight:	8 kgs	10.2 kgs
Housing dimensions:	340 x 160 x 330 mm	350 x 160 x 365 mm

(1) at -5°C

(2) at 0°C and 50 mm insertion depth

(3) e.g. 1 x for external reference probe + 1 x for probe under test (Pt 100 oder TC)

(4) optional for 115 VAC


Versions LR-Cal FLUID 100-N-2I

and LR-Cal FLUID 100-45-2I:


with 2 measuring inputs for Pt 100 and thermocouples, programmable. For unit under test and/or external reference. (Details see next page.)



Extraction of the magnetic stirrer




Calibration Bath Reservoir:



Additional at versions LR-Cal FLUID 100-N-2I and LR-Cal FLUID 100-45-2I:

Two measuring inputs for Pt 100 / Thermocouples, programmable:

- Pt 100 IEC 3-/4-wire, range -100...+700°C, accuracy $\pm 0.3^\circ\text{C}$
- Thermocouple type J, range 0...1000°C, accuracy $\pm 1^\circ\text{C}$ or types K/N/R/S, range 0...1300°C, accuracy $\pm 1^\circ\text{C}$



Included in scope of standard delivery:

- Temperature micro bath **LR-Cal FLUID 100-N** or **LR-Cal FLUID 100-45**
- Spare fuses
- Support for fixing units under test
- Connection cable for thermostat tests
- 1 Bottle (500 cm³) with silicone oil 47V10
- Cosing lids for transport purposes
- Carrying bag with shoulder strap
- Operating manual (German/English)
- Test certificate (factory certificate of calibration)



Additional at versions LR-Cal FLUID 100-N-2I and LR-Cal FLUID 100-45-2I (with 2 measuring inputs):

- Set of electrical connection cables (red/black)
- Set of clamping plugs (red/black)



Optional Accessories:


- External reference temperature sensors (see datasheet **LR-Cal LTC-F**)
- PC-Windows **software AQ2sp** incl. special RS232 connection cable. With the **AQ2sp** software, the calibrator can be completely controlled from the PC, manual or automatic calibration of one or more units under test, load of one or more test items, load and lifetime tests, creation of calibration certificates. Order-Code **590.0.000.0003.0** incl. RS232 cable.



- **Extension tube** for **LR-Cal FLUID 100-N** for increasing the immersion depth, total length 250 mm usable immersion depth 230 mm. Minimum temperature -9°C at 20°C ambient temperature. When used with silicone oil 200C5: working range $-9...+130^\circ\text{C}$, radial temperature uniformity $\pm 0.1^\circ\text{C}$ (measured 50 mm from the bottom), axial temperature uniformity $\pm 0.1^\circ\text{C}$ (or $\pm 0.15^\circ\text{C}$ at temperatures below 0°C) measured in the range 0...150 mm above the bottom. Order-Code **FLUID100-ER**.



- **Conversion** of the **LR-Cal FLUID 100-N** calibration bath into a **dry block** temperature calibrator: Temperature range: $-10...+125^\circ\text{C}$ Block made of aluminium, diameter 60 mm, useful depth 170 mm. Heating time from -10 to $+110^\circ\text{C}$: 45 min. Cooling time from 20°C to -10°C : 47 min. Stability of controlled temperature: $\pm 0.04^\circ\text{C}$. Vertical temperature uniformity: $\pm 0.03^\circ\text{C}$ at 0°C ; $\pm 0.06^\circ\text{C}$ at 80°C
 - Block without holes (for self-drilling): Order-Code **FLUID-INS-0**
 - Block with 9 holes (4.0 - 4.0 - 4.5 - 5.5 - 6.5 - 6.5 - 8.5 - 10.5 - 12.5 mm): Order-Code **FLUID-INS-9**



- **Various test liquids:** Please enquire with description of your application: dt-export@leitenberger.de, or see →



- RS232-USB interface converter